


# Bluedrop Delta Foot Scanner Quick Start Guide

Welcome and congratulations on taking the first step towards better foot health! Let's get started!



**ALWAYS KEEP OUT OF DIRECT SUNLIGHT**

1

## Position Scanner Safely on a Hard Surface

- Place scanner on the floor in your home where you can use it daily
- Keep away from sunlight and heat sources
- Do not place on thick carpet
- Leave plugged in at all times

2

## Plug Cellular Router into Power Outlet

Power LED  
should turn on



3

## Wait 30 Minutes

4

## Plug Delta Foot Scanner into Power Outlet

- Your devices will automatically connect to each other and to the internet
- The screen on the Scanner will display either of the following messages:

Step on Device

Ready

5

Note that after a few minutes the screen backlight may dim

6

## Setup is complete!

Please turn over page to learn how to scan




# Taking a Scan is Easy

Use your scanner once per day

The scanner works just like a normal weight scale


## ① Step on Scanner Barefoot

- Scanner will detect you automatically
- Position your feet as per the image on the right in the area highlighted in green.



## ② Keep still for 30 seconds

- If you move too much during your scan you might be required to take another one.



## ③ Scan Complete!

Please always leave the scanner and router plugged in to allow data to be sent



### Use your scanner once per day

- You may get a reminder text if you forget to take a scan for a couple of days



### We may contact you

- We may also check in with you periodically to ensure everything is going well