**Bluedrop Medical**

Healthcare Connected

**OneStep**<sup>™</sup>

Foot Scanner

# User Guide

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**Read instructions before use**  
**Lea las instrucciones antes**  
**de usar**

# INTENDED USE

The Bluedrop Monitoring System is intended for use by persons over the age of 18 years who are at risk of inflammation related foot conditions such as diabetic foot ulcers.

The Bluedrop Monitoring System (BMS) is comprised of the OneStep™ Foot Scanner device and its accompanying Sentinel Review Interface (SRI) software. The BMS allows the condition of the soles of the feet to be remotely monitored for thermal and visual signs of inflammation, both of which may signal the onset of certain conditions, such as diabetic foot ulcers.

Under the direction of a health care professional, the OneStep™ Foot Scanner should be used on a daily basis by the patient, in their place of residence. By placing the feet on the indicated surface from a standing or sitting position, it will measure the dermal temperature and record visual images of the soles of the feet. When used from a standing position the OneStep™ Foot Scanner will also record the patient's weight.

The SRI software will provide the data recorded by the device to a health care professional, to analyze the patient's results. The health care professional will use this data to determine if further evaluation and treatment are necessary for the patient. They will continue to monitor the images and temperature data for any persistent localized inflammation or skin changes on a regular basis and communicate with their patient as necessary.

The BMS is intended to be used under the direction of a healthcare professional as an adjunct to, and not in replacement of, self-examination and periodic foot care and examination conducted by a healthcare professional and does not diagnose any specific disease state.

## **Risks**

The OneStep™ Foot Scanner does not cause any discomfort or pain during normal use. The materials that come in contact with skin have been assessed for biocompatibility. If the device is damaged, exposed to heat or sunlight, or used improperly, inaccurate temperature readings and images may be taken which could impact your healthcare providers evaluation of potential inflammation.

## **Benefits**

Self-examination of the soles of the feet by patients for visible and thermal signs of inflammation associated with skin damage can prove difficult without assistance. Use of the Bluedrop Monitoring System (BMS) assists in the evaluation of the soles of the feet for early visible and thermal signs of inflammation in patients who are at risk of inflammation related foot conditions such as diabetic foot ulcers.

Using the data from the BMS, your healthcare provider may contact you to initiate actions to avoid potential complications like diabetic foot ulcers.

## **Description of the Device**

The OneStep™ Foot Scanner device includes sensors to collect visual and thermal information over the soles of the feet and transmit that data wirelessly to a remote server for evaluation.

The device includes a top panel on which there are an array of temperature sensors. Beneath this panel is an array of light sources and image sensors.

The user places the soles of their dry feet on the top panel of the device either by standing on the device or from a seated position. A scan is then taken by the device and is transmitted securely via the internet where it can be accessed by the user's healthcare provider for monitoring and analysis of the patient.

When the device is used from a standing position the user's weight is also recorded and displayed to the user.

## **Intended Patient Population**

The OneStep™ Foot Scanner is intended for home use by people over the age of 18 who are at risk of developing inflammatory foot diseases. The device is only suitable for use by patients who weigh 330lbs (150kg) or less. There are no specific educational requirements.

## **Contraindications**

The OneStep™ Foot Scanner should not be used by patients under the age of 18, patients weighing over 330lbs (150kg), or those who have active skin damage/ulceration on the soles of their feet.

## **Intended Operator**

The patient as recommended by the Healthcare Provider is the intended operator. A patient's care assistant may also be an intended operator. Careful handling in line with any electronic device is recommended. The patient can use all functions of the device and is not expected to carry out any servicing or maintenance except for periodic cleaning as outlined in this document.

### **Contained in the box:**

- ✓ OneStep™ Foot Scanner
- ✓ Power supply
- ✓ User guide
- ✓ Cellular router

# WARNING:

- Continue to perform your usual daily foot care routine as instructed by your healthcare provider. This will involve examining your feet for signs of developing ulcers or open wounds. Using the Bluedrop Medical Delta device is not a replacement for routine foot care.
- The maximum allowable weight is 330lbs (150kg). Use by patients over this weight can result in injury.
- Do not use this foot scanner if you are under the age of 18.
- Service of the device is not permitted. Do not carry out any modifications or maintenance of the scanner or its safe operation cannot be guaranteed and may result in injury to the user.
- Do not use with a multi socket adaptor (power strip) or an extension cord.
- Only use the included power supply (XP Power ACM18US05) to power the scanner.
- Do not use when feet or scanner are wet as the surface of the scanner may become a slip hazard.
- Ensure scanner is placed on a solid, flat, stable surface. Using the scanner on an unstable surface could result in injury.
- Do not use the scanner if you suspect it to be damaged as this may result in injury.
- Do not drop the scanner as it could cause damage to the device.
- Ensure that the power lead is not a trip hazard.
- If balance is impaired or you are at an increased risk of falls do not use the scanner from a standing position.
- The scanner may also be used from a seated position for patients who are at increased risk of falls.
- If you notice broken skin or a developing ulcer, discontinue use of the scanner and contact your healthcare provider immediately.
- Do not place the scanner in a bathroom or other humid environment, or in an area where it is likely to get wet. Moisture can degrade performance and a wet surface can become a hazard.

**330lbs/150kg Max Weight**

# Always keep out of direct sunlight

This scanner contains calibrated sensors that can be damaged by heat and sunlight

## CAUTION:

- Maximum allowable foot size is:

Region	Foot (Shoe) Size
US	14.5 (Men) 15.5 (Women)
UK	13.5 (Men and Women)
Eur	48.5 (Men and Women)

- Use in normal conditions. Do not use in temperatures outside of the range 50-95°F (10-35°C).
- Keep out of direct sunlight – sunlight can reduce sensor performance.
- Do not place scanner beside heat sources such as radiators, heaters, or vents.
- Do not place on thick carpet. Weight measurement accuracy may be reduced if placed on a carpet.
- Take care when removing scanner from packaging. Dropping the device may cause damage.
- Do not unplug the scanner while in use. Doing so may cause the device to malfunction.
- The scanner should be used for the prescribed patient only. It should not be used by more than one person.
- Do not use in combination with external cardiac defibrillator or high frequency surgical equipment or near strong magnetic field devices such as MRI. Doing so may damage the scanner.

**Bluedrop Medical**  
Healthcare Connected

# INSTRUCTIONS FOR USE

- The OneStep™ Foot Scanner allows your healthcare provider to remotely monitor your feet for early signs of skin damage.
- The scanner will take images and record the temperature of the soles of your feet over a 30 seconds period.
- For best results/monitoring it is designed to be used once per day.

# The next pages will outline some simple steps to get started...

- 
- Position the foot scanner safely in the home
  - Plug in your cellular router and wait 30 minutes
  - Plug in the foot scanner
  - Stand on the foot scanner for 30 seconds
  - You are done! You can go about your day**  
Please leave the scanner and router plugged in to allow data to be sent

Please contact technical support for assistance: (855) 453-8707


# Initial Foot Scanner Set Up

These steps only have to be done the first time you use the scanner

## Position the foot scanner safely in your home

Place the scanner on the floor in an area of your home where it is unobstructed by other objects and where you can comfortably stand or sit for at least 30 seconds per day.

It is recommended to place the scanner in the bedroom to allow it to be used when getting dressed in the morning. Do not place the scanner on thick carpet.


**Always keep  
out of direct  
sunlight**

Do not place the scanner in direct sunlight or beside a heat source. Sunlight and heat can damage the sensors and cause incorrect readings.

# Bluedrop OneStep™ Foot Scanner Quick Start Guide

Welcome and congratulations on taking the first step towards better foot health! Let's get started!

- 1 Position Scanner Safely on a Hard Surface**
  - Place scanner on the floor in your home where you can use it daily
  - Keep away from sunlight and heat sources
  - Do not place on thick carpet
  - Leave plugged in at all times
- 2 Plug Cellular Router into Power Outlet**

Power LED  
should turn on
- 3 Wait 30 Minutes**
- 4 Plug OneStep™ Foot Scanner into Power Outlet**
  - Your devices will automatically connect to each other and to the internet
  - The screen on the Scanner will display either of the following messages:

**Step on Device**

**Ready**
- 5 Note that after a few minutes the screen backlight may dim**
- 6 Setup is complete!**

Please turn to next page to learn how to scan



# Taking a Scan

## Step onto the scanner barefoot

 The maximum allowable weight is 330lbs (150kg).

Now that your scanner is set up, to use it you simply need to stand on it barefoot. You will not need to set it up every time.

- We recommend you to use the scanner every day, in the morning or as instructed by your healthcare provider
- The scanner works in the same way as a normal weight scales. You can step on the scanner when the screen is blank. The scanner will automatically detect your weight and start a scan.
- You need to ensure that each foot is placed in the sensing area (highlighted in green below). If your foot is not in this region, the scanner will not be able to take a scan of your feet.


# Taking a Scan is Easy

Use your scanner once per day

The scanner works just like a normal weight scale

## ① Keep still for 30 seconds

- If you move too much during your scan you might be required to take another one.

## ② Scan Complete!

**Please always leave the scanner and router plugged in to allow data to be sent**

### Use your scanner once per day

- You may get a reminder text if you forget to take a scan for a couple of days

### We may contact you

- We may also check in with you periodically to ensure everything is going well

- The scanner can also be used from a seated position. Place your feet in the required locations, ensuring good contact with the glass surface. Do not move your feet while the scan takes place. If the scan does not start apply a small bit more downward force with your legs.
- To ensure accurate data is collected it is important that you remain still for the full 30 seconds or until notified. If you accidentally move during the scan it may take longer, or you may be asked to remain still while the scan is repeated.

# Do not unplug the scanner or router after use

Your scanner will process and transmit data in the background and then enter sleep mode.

Processing may not take place directly after the scan is taken. During this analysis, the screen backlight will dim, and the message “Processing/Try Later” will be displayed on the screen. Once processing is complete the scanner will automatically enter sleep mode. The analysis process can take 20 minutes or longer. After the scanner has sent the data it will enter sleep mode.

## Automatic shutdown procedure.

The scanner will automatically enter sleep mode once all scan analysis and transmission is complete. We recommend leaving the scanner plugged in however If you need to unplug it, please ensure that it is in sleep mode. You can tell the scanner is in sleep mode when the screen is blank and there are no lights on. While in sleep mode the scanner can be re-activated by simply standing on it.

## Remote Monitoring

Your healthcare provider will review the data from your scanner periodically. They will be able to view photographic images of your feet, thermal information, weight data and general usage data.

This allows the healthcare provider to analyse areas of elevated temperature which can indicate inflammation. They may contact you for further evaluation based on their observations.

## Care

Once per month, gently clean the scanner surface using a slightly damp cloth and water. Dry the surface of the scanner with a dry cloth after cleaning.

## Storage

Do not store scanner in direct sunlight or near a heat source.

The scanner should be left in place between scans to allow easy daily use. We recommend leaving it in the bedroom.

The scanner should not be exposed to temperatures above 104°F (40°C) or below 41°F (5°C) during storage or to environments outside the range of 10-90% relative humidity.

# Customer Support

Phone: (855) 453-8707

Contact customer support for assistance, if needed, for:

- Setting up, using or maintaining the scanner, or
- Reporting unexpected operation or events.

## Manufacturer's Contact Details

**Website:** [www.bluedropmedical.com/support](http://www.bluedropmedical.com/support)

**Address:** Bluedrop Medical Ltd.  
Ground floor, Block 4  
Ballybrit Business Park  
Galway, Ireland H91 A4XW

# Explanation of Symbols

	Indicates the legal medical device manufacturer <b>ISO 15223-1 5.1.1</b>		Indicates the date the device was manufactured <b>ISO 15223-1 5.1.3</b>		Indicates a carrier that contains unique device identifier information <b>ISO 15223-1 5.7.10</b>
<b>IP21</b>	Mechanical ingress protection from objects >15.5mm diameter. Fluid Ingress Protection against vertical dripping water <b>IEC60529:2001</b>	 OPERATING	The highest and lowest temperature to which the product can be safely exposed during operation <b>ISO 15223-1 5.3.7</b>		Do not dispose of device in trash. Dispose the product and packaging in accordance with your local regulations. <b>BS EN 50419 Fig 1</b>
	Type BF Applied part symbol <b>ISO 7000/IEC60417 5333</b>		Serial number <b>ISO 15223-1 5.1.7</b>		Read instructions for use before operating equipment <b>ISO 7010-M002</b>
	Centre positive polarity on DC power connection <b>IEC 60417 5926</b>		Indicates the need for the user to consult the instructions for use <b>ISO 15223-1 5.4.3</b>		Not suitable for use near magnetic resonance imaging systems <b>ASTM F2503-13 Fig 9</b>
 SHIPPING	The highest and lowest temperature to which the product can be safely exposed during shipping <b>ISO 15223-1 5.3.7</b>	 STORAGE	The highest and lowest temperature to which the product can be safely exposed during storage <b>ISO 15223-1 5.3.7</b>		The highest and lowest humidity to which the product can be safely exposed during transportation and storage <b>ISO7000: 2019-2620</b>
	The power supply meets the safety requirements specified for Class II equipment <b>IEC60417: 2002-5172</b>	<b>Rx ONLY</b>	USA Federal law restricts this device to sale by or on the order of a physician <b>21 CFR 801</b>		Unique Device Identifier (UDI)
	Indicates the device needs to be protected from moisture <b>ISO 15223-1 5.3.4</b>		Keep away from heat and direct sunlight <b>ISO 15223-1 5.3.2</b>		Maximum mass/weight
	Atmospheric pressure limits <b>ISO 15223-1 5.3.9</b>		Caution <b>ISO 15223-1 5.4.4</b>		Warning <b>ISO 7010 W001</b>

Standards referenced in this document are the latest published version which is acknowledged as state if the art.

# Explanation of Common System Messages

System Message	Meaning
STEP ON DEVICE / READY	Foot scanner is ready to be used. Step on the scanner to start a scan.
ERROR / CONTACT SUPPORT	An error has occurred with the scanner. Please contact customer support via the number listed in this document or on <a href="http://www.bluedropmedical.com/support">www.bluedropmedical.com/support</a> .
NO WIFI / CHECK WIFI CONNECTION	The scanner cannot detect any Wi-Fi networks in range. Please check your home Wi-Fi network to ensure it is active or make sure you are in range.
DEVICE TOO WARM / MOVE AWAY FROM HEAT	The scanner has detected it is close to a heat source or in sunlight. Please move the scanner away from the heat source and allow its temperature to return to normal. This can take up to 2 hours.
DEVICE TOO COLD	The scanner has detected that it is below the recommended operating temperature. Please move the scanner. to a location where it is warmer and allow it to heat up to operating temperature. This can take up to 2 hours.
MAX WEIGHT / OVER 330LB / 150KG	<p>The scanner has detected that a load of over 330lbs (150kg) has been placed on the scanner.</p>  <p>WARNING: The device must not be used by people who weigh over 330lbs/150kg</p>
STAY STILL	Please remain still. A scan is in progress. If the scanner detects movement the scan process may have to restart.
SET UP WIFI / CHECK USER GUIDE	No Wi-Fi connection has been set up with this scanner. Please follow the instructions in this document to connect to your home Wi-Fi network.
SETUP COMPLETE	A Wi-Fi connection has been made successfully
SCAN COMPLETE	The scan has been completed successfully.
SCAN FAILED	The scan has failed to complete. Please try again.
SCAN CANCELLED	The user has cancelled their current scan.
PROCESSING	The scanner is currently processing data and cannot currently be used. Please try again later. This process can take 20 minutes or longer.

# Technical Information

<b>Product Name</b>	<b>OneStep™ Foot Scanner (DFS001)</b>
<b>Dimensions</b>	433.1 x 338.1 x 79.5mm <b>Net weight:</b> 5.5kg / 12lbs
<b>Max Safe Working Load</b>	150kg / 330 lbs <b>Expected service life:</b> 3 years
<b>Max Foot Size</b>	US 14.5, Eur 48.5, UK 13.5, Mondo 310
<b>IP Rating</b>	IP21 <b>Applied Part:</b> Type BF
<b>Essential Performance</b>	No essential performance
<b>MR Environment</b>	This device is not suitable for use in a magnetic resonance environment
<b>Operating Environment</b>	Temperature: 50-95°F (10-35°C) Humidity: 10-90% RH Atmospheric Pressure: 700-1060 hPA
<b>Storage Environment</b>	Temperature: 41-104°F (5-40°C) Humidity: 10-90% RH Atmospheric Pressure: 700-1060 hPA
<b>Shipping Environment</b>	Temperature: 5-104°F (-15-40°C) Humidity: 10-90% RH Atmospheric Pressure: 700-1060 hPA
<b>Power Supply</b>	Only use the accompanying power supply: XP Power ACM18US05 Input: 100-240 VAC, Frequency: 50-60Hz, Output Current: 2.5A, Output: 5V DC Class II Equipment
<b>Power (Delta Device)</b>	Input: 5.0VDC, Input Current: 2.5A

Temperature Measurement	
<b>Temperature Units</b>	°C
<b>Range</b>	15°C-35°C
<b>Accuracy</b>	±0.6°C
<b>Resolution</b>	0.1°C
<b>Precision</b>	±0.3°C

Weight Measurement	
<b>Measurement Units</b>	Kg / lbs
<b>Range</b>	30-150kg / 66-330lbs
<b>Accuracy</b>	±7.5kg / 16.5lbs
<b>Resolution</b>	0.1kg / lbs
<b>Imaging</b>	Resolution: 20-megapixel, Colour: 24bit

## Note

- Temperature sensor accuracy may be impacted if scanner is used outside its normal operating temperature range.
- Weighing function accuracy may be impacted if scanner is used outside of specified range or on a thick carpet.
- Weighing function cannot be used from a seated position.
- If the scanner has been stored at a temperature above or below the operating temperature (50 - 95°F, 10°C - 35°C), please allow at least 2 hours before use for the scanner to reach ambient temperature.

# Technical Information

## Electromagnetic Susceptibility

This device meets the requirements of the Electromagnetic Compatibility (EMC), pursuant to the Collateral Standard, IEC/EN 60601-1-2, which addresses EMC in North America, Europe and other global communities. This includes immunity to radio frequency electric fields and electrostatic discharge, in addition to the other applicable requirements of the standard. Compliance with EMC standards does not mean a device has total immunity; certain devices (cellular phones, pagers, etc.) can interrupt operation if they are used near medical equipment. Follow institutional protocol regarding the use and location of devices that could interfere with medical equipment operation.

**Note:** This device is classified as Class 1 Type BF medical electrical equipment and the device complies with specified safety levels for electrical isolation and leakage current. The XP Power AC/DC adapter has no connection to earth ground because the necessary level of protection is achieved through the use of double insulation.

### EMC WARNINGS:

- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed - to verify that they are operating normally.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the OneStep™ Foot Scanner part number, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

	Min Freq.	Max Freq.
WIFI Receiver Band	2.412GHz	2.472GHz
WIFI Transmitter Band	2.412GHz	2.472GHz

# Technical Information

Standard	Description	Electromagnetic Environment Guidance	Equipment Class
EN55011:20009 +A2:2010 ANSI C63.4:2014 EN55032:2015	Radiated Emissions	The OneStep™ Foot Scanner is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.	Group 1B
EN55011:20009 +A2:2010 ANSI C63.4:2014 EN55032:2015	Conducted Emissions	The OneStep™ Foot Scanner is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.	Group 1B
EN61000-3-2:2014 EN55032:2015	Harmonic Current Emissions	The OneStep™ Foot Scanner is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.	Class A
EN61000-3-3:2013	Voltage Fluctuations and Flicker	The OneStep™ Foot Scanner is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.	N/A

Standard	Description	Electromagnetic Environment Guidance	Test Level
EN61000-4-11:2004	Voltage Dips and Interruptions	The OneStep™ Foot Scanner (OFS) may restart as a result of power supply interruptions. If the user requires uninterrupted operation, it is recommended that the OneStep™ Foot Scanner be powered from an uninterruptible power supply.	100% dip for 10ms 100% dip for 20ms 30% dip for 500ms 100% interruption for 5s
EN61000-4-5:2014	Surge Transients	Mains power line quality should be that of a typical residential environment.	1kV (L-N) 2kV (L-E) 2kV (N-E)
EN61000-4-8:2010	Magnetic Fields	Power frequency magnetic fields should be at levels characteristic of typical location in a residential environment.	30 A/m

# Technical Information

Standard	Description	Electromagnetic Environment Guidance	Test Level
EN61000-4-3:2006 + A1:2010 Incorporating Corrigendum October 2009	Radiated Immunity	<p>Portable and mobile RF communications equipment should be used no closer to any part of the OneStep™ Foot Scanner, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended Separation Distance <math>d = [1.17] \sqrt{P}</math> <math>d = [1.17] \sqrt{P}</math>... 80MHz to 800MHz  <math>d = [2.33] \sqrt{P}</math>... 800MHz to 2.5GHz where P is the maximum output power rating of the transmitter in watts (w) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 	3V/m
EN61000-4-2:2009	Electrostatic Discharge	Mains power line quality should be that of a typical residential environment.	+/- 8kV direct contact +/- 15kV air contact
EN61000-4-4:2012	Fast Transients	Mains power line quality should be that of a typical residential environment.	+/- 2kV for AC input power lines
EN61000-4-6:2009	Conducted Immunity	<p>Portable and mobile RF communications equipment should be used no closer to any part of the OneStep™ Foot Scanner, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended Separation Distance <math>d = [1.17] \sqrt{P}</math> <math>d = [1.17] \sqrt{P}</math>... 80MHz to 800MHz  <math>d = [2.33] \sqrt{P}</math>... 800MHz to 2.5GHz where P is the maximum output power rating of the transmitter in watts (w) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 	3 Vrms 150kHz to 80Mhz 6 Vrms for Spot Frequencies

# Technical Information

## Disposal

Dispose the product and packaging according to your local regulations.

## Unauthorised Modification

No modification of this equipment is allowed.

## Maintenance

No maintenance of the scanner is required except for periodic cleaning. The user should not carry out any further maintenance of the device.

## Power Supply

The device is packaged with an external power supply (XP Power - ACM18US05). This is the only power supply that is recommended for use with this device.

## Hardcoded Wi-Fi

For testing and operational purposes, the device has been designed to automatically connect to certain Wi-Fi networks automatically using hard coded Wi-Fi details.

## Travel

The device can be used while travelling. Please note that you will need to configure the Wi-Fi in your new location.

## Effects of Pets and Pests

Keep away from pets and pests as they may inadvertently disconnect the scanner from its power supply. Scans may be taken if they stand or sit on the scanner and they exceed the minimum weight required to start a scan.

## Children

Children should not be left unsupervised in the same location as the scanner if there's a risk that they may inadvertently disconnect it from its power supply, sit on or stand on the scanner.

# Troubleshooting

# Troubleshooting

## Keep away from sources of heat - scanner too hot or cold

Do not position the scanner close to a heater, radiator, fireplace, vent or any source of heat. Do not use the scanner when it is outside of the temperature range 50-95°F (10-35°C). Exposure to heat can impact the performance of the temperature sensors. The scanner has inbuilt temperature sensors and can recognise when it is placed near a heat source. If temperature thresholds are exceeded the scanner screen will state that it needs to be moved away from the heat source or out of direct sunlight.

## Do not place in direct sunlight

**It is critical that the scanner is kept out of direct sunlight.** Sunlight can damage the scanner and can interfere with temperature sensing and visual imaging functionality. Ensure the device is placed in an area of shade to ensure it is out of sunlight throughout the day. Note the path of sunlight can change with the seasons.

## Connecting to the Wi-Fi

If you are having difficulty in connecting your scanner to the internet, you can connect your scanner to your home Wi-Fi network. To do this, input your Wi-Fi details by unplugging the scanner and starting the instructions from Step 1 in the home WiFi set-up instructions (pages 28-30). The screen may not state “Set up Wi-Fi” but you will still be able to connect to the Bluetooth Wi-Fi network like in Step 3.

## Do not unplug the device until scans have been processed

After a scan is taken it needs to be processed for 35 min prior to transmission. If the scanner is unplugged the scan will not be processed and transmitted. The scanner consumes very low power when it is operating and it is advised to leave it plugged in to ensure data transmission and processing complete successfully.

# Troubleshooting

## Wi-Fi router signal issues

- If the scanner cannot find an active Wi-Fi network, it will display the message “NO WIFI”.
- If you see this message you should check network Cellular router to make sure it is active and the power LED is on, and that the scanner is in range. Make sure the router antenna are pointed upwards.
- You can try moving the scanner closer to your Cellular router. You can move the router near to a window for a stronger signal. If you are still experiencing problems, please contact customer support.
- Keep the cellular router plugged in. If the cellular router is unplugged, unplug the scanner and begin the setup instructions on page 11 again to ensure the scanner has a connection.

## Issues with scanning

- It is important that you remain still while a scan is taking place.
- If the scanner detects movement it may have to repeat or stop the scan.
- If you are having trouble remaining still while using the scanner, you can try to use it from a seated position.
- If an error message is displayed, you should try to restart the scanner by unplugging it and re-plugging it in. If the error is still present, please contact customer support.

## Issues when using from a seated position

- Ensure as much of your foot is in contact with the glass panel as possible.
- If the scan does not start when you place your feet on the scanner, try leaning forward slightly or to press downwards with your feet. The scanner needs to feel weight before starting a scan.
- When using the scanner from a seated position it is important to note that your weight will not be displayed to you.
- You must make sure that you remain still as any movement will be detected by the scanner and may result in a scan taking longer than it should.

# TO CONNECT DIRECTLY TO HOME WI-FI

If you are having difficulty in connecting your scanner to the internet, it is also possible to connect to your home Wi-Fi using the instructions below. Please contact customer support if any assistance is required

## Step 1: Get Wi-Fi Password


Your scanner can be connected to your home Wi-fi network if you have one available.

Before we start, please have your home Wi-Fi details to hand. The password is often on the back of your Wi-Fi router.

If you have any trouble connecting, please visit [bluedropmedical.com/support](http://bluedropmedical.com/support) or contact us via our support number (855) 453-8707

## Step 2: Plug in the scanner

First, unplug both the scanner and the cellular router. Then plug only the scanner back in.

**Once you plug in the scanner you will have 10 minutes to complete the next step**

If you run out of time or wish to start again simply unplug the scanner and plug it back again to start again.

# TO CONNECT DIRECTLY TO HOME WI-FI

## Step 3: Connect to the Bluedrop Wi-Fi network

On a smartphone/laptop, connect to the Bluedrop network. You do this in the same way you would connect to any Wi-fi network.

This phone is only needed once, to connect to set up Wi-fi. You can use a phone belonging to a friend or family member if needed.


### Apple iPhone

From your home screen, go to **Settings**, then **Wi-Fi**. **Turn on Wi-Fi**. Your iPhone will automatically search for available Wi-fi networks. **Tap the Bluedrop network**.

### Android device (Samsung, Google, Motorola)

Navigate to **Settings**, then **Wi-Fi and internet** then **turn on Wi-Fi** (if its is not already on), then tap **Wi-Fi**. Then **tap the Bluedrop network**.


# TO CONNECT DIRECTLY TO HOME WI-FI

**Step 4:** You should be taken to a new webpage.

- You may alternatively be shown a pop-up notification asking you if you want to configure this Wi-Fi network.
- If you are not taken to this website please visit the web address <http://bluedropconnect.com> in your internet browser.

If something unexpected happens or you get lost, you can start the process again by un-plugging the scanner and starting from Step 1.

The image shows a smartphone screen with the Bluedrop Medical logo at the top. Below the logo, the text reads "Select your home Wi-Fi network from the list". There is a dropdown menu with three dashes "---" and a downward arrow. Below the dropdown is a blue "Refresh" button. Underneath is a "Password" label followed by a text input field. At the bottom, there is a checkbox labeled "Show Password".

**Step 5:** Select your home Wi-Fi network from the list and input password.

Once you are at this webpage you will be able to select your home Wi-Fi network and input your password and press **Connect**.

**Wi-Fi  
Connected**

If you have entered the correct password the scanner will state “Wi-Fi Connected”

If you have input an incorrect password the scanner will state “Set Up Wi-Fi”. If this happens please start this process again from Step 3.

# Troubleshooting for connecting to home Wi-Fi

## Wi-Fi signal issues

- If the scanner cannot find an active Wi-Fi network, it will display the message “NO WIFI”.
- If you see this message you should check your home Wi-Fi network to make sure it is active and that the scanner is in range.
- You can try moving the scanner closer to your Wi-Fi router. If you are still experiencing problems, you can re-input your Wi-Fi details by unplugging the scanner and starting the instructions from Step 1.
- Note: if your scanner was previously connected to another network the screen will show “Step on Device” or “Ready”, when plugged in, but you can proceed to the next step.

**Bluedrop Medical**  
Healthcare Connected

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